 IZZY

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EN

We thank you for having chosen an appliance from the IZZY range.

Safety Interlocks

Please read these instructions carefully before using the appliance and save them for future reference. Please keep the sales receipt for guarantee purposes (pls see below “Guarantee and Customer Service”). Any use which does not conform to the instructions will absolve IZZY from any liability.

- The bathroom scale is intended for indoor use. Do not use it outdoors.
- The bathroom scale is intended only for domestic use. Any other use will cancel the warranty.
- Never use the scale for any other use than indicated.
- Your scale is electronic and thus its function may be influenced by temperature variations or by its movement.
- The scale must be placed on a flat, stable, smooth, and sheltered area. Always place the scale on a sheltered area when in use as air movement and may affect its accuracy.
- Do not overload the scale by exceeding the maximum weighing capacity of 180kg/396lb as this may damage the weighing sensor.
- Do not expose the scale to any high electro-magnetic fields such as a cordless portable phone, computer, microwave etc. as this will interfere with the functioning and accuracy of the scale.
- Do not use excessive force when pressing the touch pads and never use a sharp object to operate them.
- Do not subject the scale to excessive force.
- Do not stand the scale near an oven or other direct source of heat.
- Do not operate or stand the scale in places with high humidity, where it may get wet.
- Never use highly corrosive or abstractive cleaning agents. Never use metal scouring pads for cleaning.
- Do not move the kitchen scale when in operation.

- Do not store the scale on its side.
- Do not store anything on the scale.
- Do not store the scale where you store cleaning chemicals. The vapors from some household products may affect the electronic components of your scale.
- Do not touch the kitchen scale with wet or damp hands.
- Your scale contains sensitive electronic parts. Avoid rough treatment. Do not drop, kick, or jump on it. Treat it with care to ensure the best performance.
- Always replace all batteries at the same time. Do not mix old and new batterie.
- Do not use batteries of different brand or type together.
- In case of a possible damage, do not attempt to repair the scale by yourself. Please contact one of the authorized BENRUBI service centers.
- Any repair made by an unauthorized BENRUBI Service Center CANCELS THE GUARANTEE.
- Any misuse of the device cancels the warranty.
- For your safety, this appliance conforms, to the appliance standards, regulations and directives (Low Voltage, Electromagnetic Compatibility, Environmental etc. Directives).

SAVE THESE INSTRUCTIONS

Batteries Installation & Replacement

- Make sure you connect the batteries with the correct polarity.
- If the scale is not going to be used for a long time, remove the batteries

Batteries installation:

1. Unclip the battery cover on the upper side of the scale by pushing the clip inwards and then lifting off. Place the batteries inside the battery compartment, connecting them to the appropriate polarity.
2. Clip the battery cover back into place.

Note: During use, if the battery symbol "LO" appears on the LCD screen, it means that battery power is low and the batteries should be replaced.

Batteries replacement:

1. Unclip the battery cover on the underside of the scale by pushing the clip inwards and then lifting off. Remove the old batteries and fit the replacement (3xAAA),
2. Clip the battery cover back into place

Disposal of the old batteries:

Please refer to section for the correct disposal of the product.

Language Setting / Unit of Measurement

- You may select your desired language by repeatedly pressing the button on the side of the scale.
- Your language selections are: English, Spanish, German, Greek and Croatian.
- The scale will speak the corresponding language, once the language is selected.

NOTE: Pounds units of measurement will be spoken in English. However, Spanish, German, Greek and Croatian will be spoken kilograms.

Volume Setting

- You may adjust the level of volume by the switch on the upper side of the scale.
- To turn the voice feature off, turn the switch to “-“ until it “clicks” off.

Use of Bathroom Scale

Always place the scale on a dry and flat surface prior to weighing. Do not use it on carpets.

1. Press firmly once on the scale platform to start-up the scale. The display shows “0.0” and then turns off. The scale is ready for use.
2. Place your feet on the platform. Put your feet in the same place each time and distribute your weight evenly and stand still. Your weight will automatically appear.
3. Get off the scale.
4. The scale turns off automatically after a few seconds. This protects the battery life.

Notes:

- The scale will read “Err” if it is trying to weigh above its capacity.
- The maximum weight that can be weighed is 180kg / 396lb.
- For accurate results, make the measurement at the same time (eg. 19:00-21:00 every night)
- If you do not follow the above steps, your weight measurement will not be accurate.

Tips

- Place the appliance on a stable and dry surface when not in use. Do not storage anything on the appliance in any case.
- Do not use the bath scale on the carpet.
- For accurate weight measurement, distribute your weight evenly on the appliance and stand still. Remember to place your feet on the same spot each time you are weighing yourself.

- The weight indication may be influenced by the clothes you are wearing or by the fact that you may have eaten before measurement.
- Do not weigh with wet feet. The platform may get slip.
- In case you want to observe your weight accurately for a long time period, it is important each time to weigh yourself at the same time, without having clothes on and without moving the electronic bath scale to different place.

Error Messages

	<p>This scale has a low battery indicator. It will display "LO" and the scale says "Low Battery" or readings will grow dim, or the voice quality distorts or becomes too soft to hear, the batteries need replacement.</p>
	<p>"O-Ld" indicates an overload warning. The maximum capacity of this scale is 396 lbs / 180 kgs. Should "O-Ld" appear on the display during a weighing, this indicates that the maximum weighing capacity of the scale has been exceeded. Remove the weight immediately; otherwise, permanent damage may occur.</p>
	<p>"Err" indicates the scale is unable to get a reading. The "Err" will display along with the voice message "Error". Step off the scale and repeat the weight measurement. Your weight cannot be locked-in if you do not stand still.</p>

Care and Cleaning

- Do not immerse the electronic bath scale into water or any other liquid. Do not wet it with water or any other liquid.
- The glass platform may be pretty slippery when wet, thus please make sure both the glass your feet are dry before using.
- Please use soft tissue with alcohol or glass cleaner to clean the surface if it is dirty. No soap or other chemicals. Keep it away from water, heat and extreme coldness.
- Always protect the appliance from high temperatures and/or strong sunshine.
- Never use chemical or abrasive cleaners for cleaning the appliance.
- Handle your scale with care; do not drop it or jump on it.
- Do not step on the scale with shoes. This may cause damage of the glass platform.

Technical Characteristics

Model: Bathroom scale Cloud Talking IZ-7011

Max weight: 180kg/396lb

Graduation: 0.1kg / 0.2lb

Weight measuring units: Kilograms (Kg), Libres (Lb)

Batteries: 3xAAA (1.5V each) – included

Important Information for Correct Disposal of the Product in Accordance with EC Directive 2002/96/EC

At the end of its working life, the product must not be disposed of as urban waste. It must be taken to a special local authority differentiated waste collection centre or to a dealer providing this service. Disposing of a household appliance separately avoids possible negative consequences for the environment and health deriving from inappropriate disposal and enables the constituent materials to be recovered to obtain significant savings in energy and resources.

As a reminder of the need to dispose of household appliances separately, the product is marked with a crossed-out wheeled dustbin.

Guarantee & Customer Service

- BENRUBI Company guarantees its products for any manufacturing defect or defective materials for a period of two (2) years from the date of purchase or delivery by presenting purchase receipt. If your product malfunctions or you find any defects, please revert to the retail location you have purchased it from or contact an authorized BENRUBI Service Centre. To find the nearest authorized BENRUBI Service Centre, visit our website www.benrubi.gr. It is in the discretion of BENRUBI Company to replace or repair the defective product.

- The warranty does not cover damages resulting from normal wear, breakage, improper installation or maintenance of the product, mishandling, abnormal working conditions, failure to implement the instructions, conversion or repair of the product by an unauthorized technician who does not belong to the BENRUBI Service Centers.

Also, the warranty does not cover the following indicative factors:

- Stains, discoloration or scratches.
- Burns due to exposure to heat or flame.
- Damage from thermal shock (sudden temperature change or change of voltage).

FOR MORE INFORMATION PLEASE CONTACT:

(NEW) H. BENRUBI & FILS SA

27 Aghiou Thoma str, GR15124, Maroussi – Athens,

Tel./Customer Support: +30210 6156400

e-mail: benrubi-sda@benrubi.gr